

Covenant Health
Accu-Chek Inform II Competency Sheet

NAME:			UNIT: CSON		DATE:	
			STUDENT ID:			
INITIAL TRAINING	6 MONTH	ANNUAL ☐				
Assessment	PASS	FAIL	METHOD			
HealthStream Interactive learning tool completed			certificate attached			
Assemble testing supplies:						
Meter, strips, tote, QC solutions- high/low,			observed			
alcohol wipes, 2X2s, gloves			observed			
Applies appropriate PPE - gloves			observed			
Patient test/Simulated patient test						
Identifies patient using 2 patient identifiers			observed			
Explains procedure to patient			observed			
Turns meter on, Scans operator ID			observed			
Verbalizes/carries out process if ID fails			observed			
Selects appropriate function (patient test)			obs/tchbk			
Scans patient armband			observed			
Verifies test strip code matches vial			observed			
Inserts strip, recaps vial immediately			observed			
Correct site care and specimen aquired			observed			
Applies sample to front of strip correctly			observed			
Verbalizes/demonstrates response to abnormals/criticals (comments/ repeats)			obs/tchbk			
Confirms result/comment(s)			observed			
Cleans, performs, documents maintenance			observed			
*must stay "wet" for 1 minute			observed			
*not in strip port or on metal contacts			observed			
Can verbalize approved cleaning agents			tchbk			
Quality Control						
Verbalizes when QC must be performed			obs/tchbk			
Verbalizes LOCKOUT			obs/tchbk			
Verbalizes/Demonstrates high/low QC			obs/tchbk			
Verbalizes/demonstrates special docking station/use (and location if available)			obs/tchbk			
Results Review						
All results			observed			
Patient-specific results			observed			
Quiz Score:			Pass / Fail			
This CSON student has satisfactorily completed the competency for the Accu-Chek Inform II bedside blood glucose monitor.						
Student Signature: _____						
Validator Signature: _____						
Validator's Printed Name: Cherish Walker, MSN, RN						
This CSON student has NOT satisfactorily demonstrated competency. This form is being forwarded to _____ for further education/training. Date: _____						
Student Signature: _____						
Validator Signature: _____						
Validator's Printed Name: _____						

1. Personal protection equipment must be worn when performing blood glucose or glucose control testing and when cleaning and disinfecting the Accu-Chek Inform II meter.

True False

2. The proper site for the finger puncture is on the side of the finger.

True False

3. The Accu-Chek Inform II meter should be docked in the base unit when not in use.

True False

4. Glucose control solutions do not have to be dated when first opened.

True False

5. Testing will start when the appropriate amount of sample is applied to the Accu-Chek Inform II test strip, about 0.6 ul.

True False

6. Control solutions are stable for _____ months after opening?

a. 2 months 4 months

no expiration date 3 months

7. What type of blood cannot be used with the Accu-Chek Inform II meter and the Accu-Chek Inform II test strips?

a. Capillary Venous/Arterial

Neonatal Cord blood

8. You _____ scan the patient's armband instead of the patient's chart label to avoid errors.

a. should should not