





Expanding your clinical experience

Philips ultrasound women's healthcare
guide to sonography of the breast

PHILIPS

Guide to Sonography of the Breast



Ultrasound of the breast has been performed for more than a decade, with exam protocol continuously evolving based on image quality of the current equipment. Sonography provides a noninvasive, tomographic display of the breast without ionizing radiation.

Continuing advances in digital technology bring added benefits to ultrasound as a useful complement to mammography and physical examination in the evaluation of breast disease. The adjunctive use of high-frequency imaging ultrasound with mammography provides increased confidence in differentiating solid lesions to determine those that are benign, reducing the need to biopsy many lesions considered indeterminate on mammograms.

Breast Anatomy

Sonography of the breast requires that the sonologist or sonographer has a comprehensive knowledge of the anatomy of the breast. The sonographic presentation of the breast, which is composed of fat, fibrous tissue and glandular tissue, depends greatly on the hormonal status of the patient.

Sonographic anatomy of the breast


- Skin
- Subcutaneous fat
- Breast parenchyma (mammary zone)
- Copper's ligaments
- Retromammary fat
- Pectoralis muscle
- Ribs/pleura
- Nipple region
- Tail of Spence

This image demonstrates most of the anatomical elements encountered by breast sonography.

The skin is seen as a highly reflective band along the surface of the breast. Normal thickness is 2 mm. Subcutaneous fat lies between the skin and the parenchymal (mammary zone) tissue. The quantity of fat varies. Fibroglandular tissue is the echogenic layer of tissue beneath the nipple and subcutaneous fat. Once again the ratio of fat versus fibroglandular tissue varies among patients and may also depend on the age, parity and hormonal status of the patient.

Cooper's ligaments are suspensory ligaments, providing support to the glandular structures of the breast. They extend radially from the deep fascial planes to the skin. Retromammary fat forms a layer between the deep fascial plane of the breast and the pectoral muscle, defining the posterior boundary of the glandular tissue. The pectoralis muscles can be clearly imaged sonographically in the direction of their fibers. They appear above the ribs and parallel to the skin. The clinician must have a clear understanding of normal appearance to evaluate the area for tumor extension.

The ribs are readily identified laterally because bone attenuates the sound beam, resulting in acoustic shadowing. Medically, the ribs appear as hypoechoic structures containing low-level echoes.



The nipple region requires special attention because it consists of both dense connective tissue (of the nipple) and partially connective tissue of the lactiferous ducts, which can cause posterior acoustic shadowing.

The tail of Spence (axillary tail) is the portion of the breast that extends into the axillary region.

Technique

The examination of the breast with ultrasound requires a preview of the mammogram as well as a good physical examination and patient history.

The patient history should include:

- Family history
- Previous mammograms
- Masses
- Scars
- Skin changes
- Nipple discharge
- Breast contour
- Age
- Parity
- Gravida
- Aborta
- Medications (hormones)
- Surgeries of the breast

Upon completing a thorough review of the patient's history, mammograms, and physical examination, the breast ultrasound exam can begin. The patient is scanned supine with the ipsilateral hand either above the head, or on the hip with the elbow pointed back.

This causes the breast to flatten across the pectoralis muscle. Depending on the size of the breast, multiple scanning positions may be required. A cushion is placed behind the shoulder of the breast being examined. These techniques help stabilize the breast and provide reproducible positioning if open surgery is required.

The breast is then scanned with a high-frequency (at least 7-15 MHz) transducer. The breast may be scanned longitudinally or transversely. Stavros et al. recently described scanning of the breast in the radial and antiradial presentation. Regardless of which method is used, the breast needs to be examined in two or more orthogonal imaging planes when a mass is identified.

The breast also needs to be clearly labeled while scanning is being performed. While the "face of the clock" is the most common labeling method, others prefer to divide the breast into four quadrants. Regardless of the method used, consistency must be maintained to allow reproducibility in follow-up scans and quality control within imaging.

It is important to pay close attention to the nipple area when scanning, since a shadow can be caused by the erection of the nipple. In this case, either apply more pressure, use more gel, or scan the nipple area obliquely. Scanning the nipple region obliquely is easily performed by placing a rolled-up towel between the patient's breasts. By having her roll onto her side, the nipple of interest will be on top. The nipple can then be scanned from the side. This reduces the shadowing, avoids the need for a standoff pad, and provides adequate visualization of the anatomy posterior to the nipple.

Indications for Breast Ultrasound

Breast ultrasound is used as an adjunct to mammography and physical examination. The most common indications to perform an ultrasound exam are the presence of a palpable mass or discovery of a mass on mammogram. Ultrasound assists in identifying the mass as cystic or solid. High-frequency imaging assists in differentiating solid masses and identifying those lesions that are more likely benign and for whom a tissue diagnosis biopsy is more optional. This information enables the physician and patient to make a decision as to how to manage the mass.

Guide to Sonography of the Breast_{continued}

Ultrasound guidance of aspiration, fine needle aspiration (FNA) and core biopsies is a rapidly growing application for ultrasound of the breast. With advances in technology, ultrasound-guided needle biopsies offer both the patient and the sonologist a simple, effective choice in breast management.


Indications

- After abnormal mammogram for differentiation between cysts and solid masses
- For differentiation of solid masses to determine more likely benign lesions
- Palpable mass not visible in a radiographically dense breast
- Young, pregnant or lactating patient with a palpable mass
- Suspected abscess in infected breast
- Mass that cannot be completely evaluated with mammography because of location
- Guidance for interventional procedures

Other potential circumstances for ultrasound examination include suspected leaks from silicone implants and oncology follow-up.

Sonography is advantageous because it provides a painless, noninvasive tomographic study of the breast. The strengths of sonography complement mammography, providing additional anatomical information and increasing diagnostic confidence.

The evaluation of masses in the breast has been traditionally one of investigating cystic versus solid. Using high-frequency ultrasound, even more information is obtainable and solid masses can be evaluated based on their ultrasound characteristics. The ability to differentiate benign masses from other suspicious tissue may reduce the need for biopsy for many patients.



The breast can be scanned longitudinally or transversely whether the entire breast is scanned or just an area of interest (see Figure A). If a mass is discovered, the radial/antiradial approach can help further identify borders and the possibility of ductal extension (see Figure B).

Ultrasound Characteristics

- Margins
- Shape
- Echogenicity
- Echotexture
- Orientation
- Posterior acoustic attenuation pattern

Overview of differentiating solid breast masses with high-frequency ultrasound

Extensive use of ultrasound for adjunctive breast exams has shown that lesions have definite image characteristics that indicate benign appearance. When a solid breast mass is encountered, the following criteria should be evaluated:

- Margins: degree of irregularity
 - Benign masses usually are indicated by smooth margins.
 - Malignant tumors appear aggressive and may have finger-like extensions or spiculations.
- Shape: ovoid, irregular, lobulated or spherical
 - Benign masses usually are spherical or ovoid – a smooth, round or egg shape. Stavros has found that lesions with three or less gentle lobulations are usually benign.
 - Malignant masses tend to be variable or irregular in shape.

- Echogenicity and echotexture: characterization of echo pattern and texture
 - Benign masses usually are homogeneous, or of equal or lower echogenicity to the surrounding tissue, and have uniform internal echoes.
 - Malignancies are most often heterogeneous with a variety of echotextures, predominately hypoechoic.
- Orientation: position to tissue planes
 - Benign lesions usually grow parallel to tissue planes.
 - Malignant masses tend to invade surrounding tissue and cross tissue planes.
- Posterior acoustic attenuation pattern: degree of sound beam absorption
 - Benign solid masses have minimal through transmission and cystic lesions have high transmission of sound.
 - Malignant masses may absorb the sound beams and cause echoes posterior to their location to be very black and chaotic. Shadowing is also chaotic or random. (Colloid cancers may have through transmission.)



Simple cysts

- Completely anechoic
- Well-circumscribed
- Thinly-encapsulated
- Enhanced through transmission
- Thinly-edged shadows

Solid masses – probably benign

- Oval shape
- Smooth defined borders
- Uniformly low or medium-level internal echoes
- Minimal attenuation, if any
- Wider than tall
- Two to three gentle lobulations

Malignant masses

- Variable shapes
- Irregular, ill-defined borders
- Low-level, non-uniform internal echoes
- Taller than wide
- Posterior attenuation



CLINICAL SOURCE
Amy M Lex, MS, RT(R), RDMS, Philips Ultrasound, Bothell, WA
Paula Gordon, MD

REFERENCES

- Stavros AT, Thickman D, Rapp C, Dennis M, Parker SH, Sisney GA. "Solid breast nodules: use of sonography to distinguish between benign and malignant lesions." *Radiology* 1995;196:123-134.
- Jackson VP. "The role of US in breast imaging." *Radiology* 1990;177:305-311.
- Carr-Hoefer C: Breast. In *Diagnostic Medical Sonography: A guide to clinical practice*, Kawamura D (ed.), J.B. Lippincot: 1992:541-578.
- Fornage BD. "Ultrasound of the breast." *Ultrasound Quarterly* Vol. 11. No. 1, 1993, pp. 1-39.
- Leucht W. *Teaching Atlas of Breast Ultrasound*. Thieme Medical Publishing, Inc.: 1992.

**Philips Medical Systems is part of
Royal Philips Electronics**

Interested?

Would you like to know more about our imaginative products? Please do not hesitate to contact us.

We would be glad to hear from you.

On the web

www.medical.philips.com

Via email

medical@philips.com

By fax

+31 40 27 64 887

By mail

Philips Medical Systems
Global Information Center
P.O. Box 1286
5602 BG Eindhoven
The Netherlands

By phone

Asia

Tel: +852 2821 5888

Europe, Middle East, Africa

Tel: +49 7031 463 2254

Latin America

Tel: +1 954 628 1000

North America

Tel: +800 285 5585



© 2006 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Medical Systems Nederland B.V. reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in USA
4522 962 09961/795 * NOV 2006