

Experimental 4: Heat of Neutralization

A. Heat Capacity of Calorimeter

Assuming: heat capacity of calorimeter is zero !!!

B. Heat of Neutralization of HCl+ NaOH

	1.0 M NaOH, 50 mL	1.0 M HCl, 50 mL
Temperature before mixing	22°C	22°C
Average (T _i)	22°C	
Temperature after reaction (T _f)	32°C	
Calculate Heat of Neutralization		

C. Heat of Neutralization of CH₃COOH+NaOH

	1.0 M NaOH, 50 mL	1.0 M CH ₃ COOH, 50 mL
Temperature before mixing	22 °C	21.5 °C
Average (T _i)	22 °C	
Temperature after reaction (T _f)	27°C	
Calculate Heat of Neutralization		

D. Compare part B and C, discuss the results

Group members:

Instructor: