

SEEK & FIND: OSTOMY



Student Name Hyun You Date 1/19/21

Instructor Signature _____ Student's earned points _____/20

Please use the product name & NOT the product number when completing this assignment.	Possible Points	Earned Points
<p>Name 2 pouching systems in the listed category. A pouching system is a barrier and a pouch. Describe abdominal characteristics, stoma characteristics, and one other reason why you would chose that system.</p> <p>Include manufacturer and full product name; it should be available in your home country.</p> <p>For each category, be sure to list products from <u>two separate companies</u>.</p> <p>Example: Cymed one piece, drainable pouch with thin hydrocolloid washer and press 'n seal closure (cannot use this one as an answer).</p> <p>1. Flat, one-piece pouching system</p> <ul style="list-style-type: none"> a. Hollister Premier One-piece Drainable Pouch w/ FLAT care plus barrier b. Coloolast Sensura One-Piece Mio Flat MAXI drainable pouch <p>Abdominal and stoma characteristics/use</p> <p>Abdomen is firm and stoma is flush and/or budded. One-piece system is better for the Pt who needs simplicity or who lacks the dexterity or hand strength to snap the pouch onto the wafer.</p> <p>2. Flat, two-piece pouching system</p> <ul style="list-style-type: none"> a. Hollister Two-piece New image drainable pouch w/ new image Flat flexextend skin barrier b. Coloplast Two-piece Assura MAXI drainable pouch, opaque w/ Assura Flat skin barrier cut-to-fit. <p>Abdominal and stoma characteristics/use</p> <p>Abdoemn is firm, and stoma is flush and/or budded. Two-piece system is better for the Pt who wants to rinse and swap pouches or who has limited visual acuity.</p>	3	
<p>Name 2 pouching systems in the listed category. A pouching system is a barrier and a pouch. Describe abdominal characteristics, stoma characteristics, and one other reason why you would chose that system</p> <p>Include manufacturer and full product name; it should be available in your home country.</p> <p>For each category, be sure to list products from <u>two separate companies</u>.</p>	3	

**** This assignment may also be completed via an internet search.
2018**

SEEK & FIND: OSTOMY



<p>Example: Genairex 10" One-piece extended pre-cut wear urostomy system with built-in convexity (cannot use this one as an answer).</p> <p>1. Soft/flexible convexity</p> <ul style="list-style-type: none"> a. Hollister Two-piece New image extended convex skin barrier w/ new image drainable pouch. b. Coloplast Two-piece Sensura Xpro convex light MAXI urostomy pouch transparent with extended wear cut-to-fit barrier. <p>Abdominal and stoma characteristics/use</p> <p>Abdomen is firm and stoma is retracted, or Abdomen is soft, but stoma is budded. Soft/flexible convexity increased the wear time of the ostomy pouching system. If the Pt only afford few numbers of ostomy system, it is a good choice to use.</p> <p>3. Hard/firm convexity</p> <ul style="list-style-type: none"> a. Coloplast One-piece Assura Deep convex MAXI drainable pouch w/ cut-to-fit standard wear barrier w/ easi close wide outlet. b. Convatec One-piece Active Life Deep convex pre-cut urostomy pouch with durahesive skin barrier. <p>Abdominal and stoma characteristics/use</p> <p>Abdomen is very soft and stoma is flush and/or retracted. Obese Pts w/ liquid or loose stool output need hard convexity pouching system.</p>		
<p>Provide manufacturers' name(s) and give use for the following accessory products that is: Why would one use the accessory product?</p> <p>1. Skin barrier rings (list two)</p> <ul style="list-style-type: none"> a. Coloplast Brava Moldable Rings b. Hollister Adapt Barrier Rings <p>Use: To fill in or calulk uneven skin contours near stoma to create a flatter surface and help prevent ostomy drainage from getting under the ostomy barrier.</p> <p>2. Disposable convex insert (NOT convex barrier ring)</p> <p>Convatec SUR-FIT Natura dispoable convex insert</p> <p>Use: Disposable insert ideal for people with stoma complications such as a recessed stoma, retracted stoma, or flush stoma.</p> <p>3. Liquid bonding cement</p> <p>Torbot Bonding Cement</p> <p>Use: Covieniently adhere ostomy appliance to the skin.</p> <p>4. Ostomy Belt (NOT hernia belt)</p> <p>Coloplast Brava Ostomy Belt</p> <p>Use: Support the weight of the ostomy appliances which attach to the stoma.</p>	3	

**** This assignment may also be completed via an internet search.
2018**

SEEK & FIND: OSTOMY



<p>5. Skin barrier powder: Pathenon K-Gem powder</p> <p>Use: Heal red, excoriated skin around the Pt's stoma.</p> <hr/> <p>6. Pouch deodorant: Safe n' simple Ostomy pouch deodorant</p> <p>Use: Eliminate odor in an ostomy appliance and helps in emptying pouch content.</p>		
<p>Name 2 different pediatric pouches. Provide an alternate use for a small pouch.</p> <p>Include manufacturer and full product name; be sure to list products from <u>two separate companies</u>.</p> <p>Example: Montreal Ostomy (BE/MO) One-piece pediatric urostomy pouch (Cannot use as an answer).</p> <ol style="list-style-type: none"> 1. ConvaTec Little One-Piece Cut-to-fit Transparen drainable pouch w/ stomahesive skin barrier. 2. Coloplast Sensura Two-piece Mio Baby barrier w/ Sensura Mio Baby flex drainable pouch <p>Alternate use: For intimate moments, special pouches are available that are designed to be small and less bulky.</p>	2	
<p>Name a commercial access port used to bring a tube out through the front of a pouch. Include both manufacturer and full product name.</p> <p>Universal Catheter Access Port (UCAP) by hollister</p> <p>Your facility does not carry any commercial access ports, and you are asked to bring a tube out through a pouch. Describe an alternate method for bringing the tube out through a pouch. (The nipple device described in the textbook in the tube chapter is a commercial device; it is <u>not</u> an alternate method).</p> <p>A foley bag may be connected to the system for dependent drainage. A tube can be run through nipple device in the wound manager and connected to low continuous suction.</p>	2	
<p>Select <u>one</u> major manufacturer of ostomy supplies. List 3 types of skin barrier wafers that have different ingredients integrated into the wafer. Include an indication for use for each type of <u>skin barrier wafer</u> chosen.</p> <ol style="list-style-type: none"> 1. Hollister New image Forma Flex moldable skin barrier <p>Use: The Forma Flex barriers is an extended wear skin barrier w/ a specialized formulation that is designed to be stretched and shaped by hand.</p>	3	

**** This assignment may also be completed via an internet search.
2018**

SEEK & FIND: OSTOMY



<p>It is the possible choice when the stoma is not round or still changing size. The Pt doesn't need to carry scissor for air travel.</p> <p>2. Hollister New image Flat ceraplus skin barrier.</p> <p>Use: The CeraPlus skin barrier is infused with ceramide, a naturally occurring lipid in the skin. Ceramide keeps skin protected and maintains moisture barrier. It helps peristomal skin intact and healthy. Therefore, ostomy appliances stay in place better.</p> <p>3. Hollistr New image Flat flexWear skin barrier.</p> <p>Use: This skin barrier with floating flange protects skin from corrosive drainage. It eliminates pressure on tender post-operative abdomens with the exclusive floating flange. The Pt assures the floating flange which provides utmost security with the patented locking system.</p>		
<p>Name 2 pouches/pouching systems from different manufacturers that are specifically designed for wound/fistula needs.</p> <p>Include manufacturer and full product name.</p> <ol style="list-style-type: none"> 1. ConvaTec Eakin Fistula Wound Pouch 2. Coloplast Fistula and wound management system 	2	
<p>List 2 pouching systems for large stomas that have skin barriers with <i>greater than a 4 inch (10.16 cm) cutting surface</i>. Do <u>not</u> include a fistula pouch (see prior question). Include manufacturer and full product name.</p> <p>Example: Torbot Bongort® Max E Wound Management Pouch with a ConvaTec Stomahesive skin barrier 8 x 8 (cannot use as an answer).</p> <ol style="list-style-type: none"> 1. Hollister Wound Drainage Collector w/ Hollister Flexend skin barrier 8x8 2. ConvaTec Eakin Wound Pouch w/ Colopanplast skin barrier 8x8 	2	

**Submit completed project to your dropbox.

9/2020

**** This assignment may also be completed via an internet search.
2018**