

**SEEK & FIND: OSTOMY**



Student Name: Stacy Mariano\_\_\_\_\_ Date 4/7/2023\_\_\_\_\_

Instructor Signature \_\_\_\_\_ Student's earned points \_\_\_\_\_/20

This assignment is intended to help you become familiar with some of the products used when providing ostomy care. Carefully read the information to help you identify the products that meet the descriptions being provided.

\*\*It is important to note a pouching system is a skin barrier wafer and a pouch. A complete answer should include both, unless otherwise indicated. Include the manufacturer and full, product name. Product numbers should not be used.

	Possible Points	Earned Points
<p>Name 2 pouching systems in each of the listed categories. The products should be from two separate companies.</p> <p>Example: Cymed one-piece, drainable pouch with thin hydrocolloid skin barrier (washer) and press 'n seal closure. This example can not be used as one of your answers.</p> <p>1. Flat, one-piece pouching system</p> <p>a. Convatec 1 pc Drainable pouch w/moldable technology (transparent or opaque).</p> <p>b. Coloplast Mio Extended Wear easy close, CTF.</p> <p>_____</p> <p>Describe abdominal characteristics, stoma characteristics, and one other reason why you would choose this system. For abdomen that is mostly flat with firm supportive tissue and budded stoma. Also would choose for ease of application.</p> <p>_____</p> <p>2. Flat, two-piece pouching system</p> <p>a. _Hollister New Image Flexwear Standard Flat.</p> <p>b. _Coloplast Sensura Mio Flex.</p> <p>**Both products come in pre-cut or CTF.</p> <p>Describe abdominal characteristics, stoma characteristics, and one other reason why you would choose this system. For abdomen that is flat with budded stoma, CTF for varied stoma shape, good for patients who change pouch frequently or have excessive gas.</p> <p>_____</p>	3	

**SEEK & FIND: OSTOMY**



<p>Name 2 pouching systems in each of the listed categories. The products should be from two separate companies.</p> <p>Example: Genairex 10" one-piece extended wear, pre-cut urostomy system with built-in convexity. This example can not be used as an answer.</p> <p>1. Soft/flexible convexity</p> <p>a. <u>Marlen 1 piece Convex with Aquatack shallow. Pre-cut convexity.</u></p> <p>b. <u>Coloplast Sensura Mio Flex, Convex light, CTF or pre-cut.</u></p> <p>Describe abdominal characteristics, stoma characteristics, and one other reason why you would choose this system.</p> <p><u>Use for semi-firm support structures, low profile or flat stoma. May also be used for patients with loose surrounding skin.</u></p> <p>3. Hard/firm convexity</p> <p>a. <u>Convatec Surfite Natura Durahesive Pre-cut barriers with Convex-It.</u></p> <p>b. <u>Coloplast Sensure Mio Flex Convex Deep 1 piece. CTF</u></p> <p>Describe abdominal characteristics, stoma characteristics, and one other reason why you would choose this system.</p> <p><u>Use with soft abdomen and supportive tissues, flush or retracted stomas, and deep crevasses. May also be used for patients with an abundance of scarring on the abdomen. The Convatec product is good for liquid output and the Coloplast Mio boasts superior flexibility for active lifestyles.</u></p> <p>_____</p> <p>_____</p>	3	
<p>Below is a list of accessory products. For each category, identify products from different manufacturers. Identify <u>why</u> the accessory products would be used.</p> <p>1. Skin barrier rings</p> <p>a. <u>Coloplast Brava moldable rings.</u></p> <p>b. <u>Hollister Adapt CeraRing Oval Convex Barrier ring.</u></p> <p>Use: to fill in gaps, provide convexity, adaptable for irregular shaped stomas.</p> <p>_____</p> <p>2. Disposable convex insert (NOT convex barrier ring)</p> <p><u>Convatec SurFit Natura Convex Insert</u></p> <p>Use: light convexity for firm skin or slightly budded stoma. _____</p> <p>3. Liquid bonding cement <u>Torbot Liquid Bonding Cement. (latex, flammable)</u></p>	3	

**SEEK & FIND: OSTOMY**



<p>Use: Extra strong adhesive for patients with problems with adherence.</p> <hr/> <p>4. Ostomy Belt (NOT hernia belt) _Hollister Ostomy belt (2 sizes).</p> <hr/> <p>Use: Extra security especially for patients with abdominal crease at waistline when bending.</p> <hr/> <p>5. Skin barrier powder _Coloplast Brava stoma powder.</p> <hr/> <p>Use: moisture absorption, crusting for denuded skin protection.</p> <hr/> <p>6. Pouch deodorant _Hollister Adapt deodorant/lubricant _____</p> <p>Use: _neutralize odors and helps effluent move from stoma to pouch, avoiding backup. _____</p>		
<p>Name 2 pediatric pouching systems from <u>two separate companies</u>. Include manufacturer and full product name.</p> <p>Example: Coloplast Assura AC pediatric one piece drainable pouch. This example can not be used as one of your answers.</p> <p>1. _Hollister Pouchkins 2 piece drainable pouch.</p> <hr/> <p>2. _Convatec Little Ones 2 piece drainable pouch system.</p> <hr/> <p>Name one alternative use, other than for a pediatric individual, for a pediatric pouch. _May use for someone with ultra-sensitive skin.</p> <hr/>	2	
<p>Access ports help to bring tubes through an ostomy pouch. Name a commercial access port. Include both manufacturer and full product name.</p> <hr/> <p>When an a commercial access port is not available, how might one bring a tube through a pouch and secure the tube? Describe an alternate method for bringing the tube out through a pouch <u>without</u> using an access port or a commercial nipple device. Be sure to describe the procedure in a manner others can replicate the process.</p> <hr/> <hr/>	2	

**SEEK & FIND: OSTOMY**



<p>Select <u>one</u> major manufacturer of ostomy supplies (Coloplast, Convatec, Hollister, NuHope, etc). List 3 types of <u>skin barrier wafers</u> that have <u>different ingredients integrated into the wafer</u>. Identify when each skin barrier wafer would be chosen for use. For example, hydrocolloid skin barrier wafer.</p> <p>1. Convatec: Sure Fit Duraheasive wafer. _____</p> <p>Use: absorbs moisture and swells up to create secure seal. Use with liquid effluent. _____</p> <p>2. Hollister Ceraplus wafer contains ceramide to eliminate moisture loss from damaged skin. _____</p> <p>Use : Use with macerated or denuded skin, delicate skin, or sensitive skin. Latex free _____</p> <p>3. _Coloplast Sensura 2 layer technology with strong adhesive to provide secure seal and skin layer to protect skin by absorbing moisture. _____</p> <p>Use __Flexible barrier wafer for active lifestyles. _____</p>	3	
<p>Name 2 pouching systems from different manufacturers that are specifically designed for wound/fistula needs. Include manufacturer and full product name.</p> <p>1. Convatec Eakin Fistula Wound Pouch with New Tap Closure.</p> <p>2. Coloplast Fistula Wound Management System.</p>	2	
<p>Identify 2 pouching systems with a skin barrier wafer that allows for cutting an opening <u>larger than 4 inches (101.6 cm)</u>. Do <u>not</u> include a fistula pouch (see prior question). Include the name of the manufacturer of the pouching system and full product name.</p> <p>1. _Coloplast Sensura PostOp Ostomy, 1 piece. _____</p> <p>2. Hollister 12 inch Premier CTF, Hi Output Ultra Clear, 1 piece. _____</p>	2	
<b>Total Points</b>	20	

**Submit completed project to your dropbox.**