

**Strip 5-1**

Rhythm: Regular  
Rate: 79  
P waves: Sinus  
PR interval: 0.14 to 0.16 seconds  
QRS: 0.06 to 0.08 seconds  
Comment: An inverted T wave is present

**Strip 5-2**

Rhythm: Regular  
Rate: 45  
P waves: Sinus  
PR interval: 0.14 to 0.16 seconds  
QRS: 0.08 seconds

**Strip 5-3**

Rhythm: Regular  
Rate: 88  
P waves: Sinus  
PR interval: 0.20 seconds  
QRS: 0.08–0.10 seconds  
Comment: A depressed ST segment and biphasic T wave is present

**Strip 5-4**

Rhythm: Irregular  
Rate: 50  
P waves: Sinus  
PR interval: 0.16 to 0.18 seconds  
QRS: 0.04 seconds

**Strip 5-5**

Rhythm: Regular  
Rate: 50  
P waves: Sinus  
PR interval: 0.16–0.18 seconds  
QRS: 0.06–0.08 seconds  
Comment: An elevated ST segment is present

**Strip 5-6**

Rhythm: Regular  
Rate: 136  
P waves: Sinus  
PR interval: 0.14 to 0.16 seconds  
QRS: 0.06 to 0.08 seconds

**Strip 5-7**

Rhythm: Regular  
Rate: 68  
P waves: Sinus  
PR interval: 0.16 to 0.18 seconds  
QRS: 0.12 to 0.14 seconds  
Comment: A U wave is present

**Strip 5-8**

Rhythm: Irregular  
Rate: 50  
P waves: Sinus  
PR interval: 0.12 to 0.14 seconds  
QRS: 0.06–0.08 seconds  
Comment: An elevated ST segment and inverted T wave are present

**Strip 5-9**

Rhythm: Regular  
Rate: 94  
P waves: Sinus  
PR interval: 0.14 to 0.16 seconds  
QRS: 0.06 to 0.08 seconds  
Comment: A depressed ST segment is present